#### Database surface electromyography (SEMG) for lower limb analysis

**Abstract:** This database contains samples from 11 subjects with knee Abnormality previously diagnosed by a professional and 11 normally. These data were collected with electromyography and goniometry equipment MWX8 Datalog Biometrics.

Característica	Time Series	Number of	22	Area	Señales físicas
de la base de		samples			
datos					
Atributes	Real	Number of	5	Sampling	July 2012
		attributes		date	March 2013
		attributes		dute	July 2013

#### 1. Sources:

Batallón de sanidad (BASAN) con apoyo de la Universidad Militar Nueva Granada – Bogotá (samples July 2012).

Carrera 7 No 52-48, Bogotá.

TecnoParque SENA nodo Manizales (samples march - july 2013).

Km 10 route to Magdalena, Manizales.

#### 2. Information database:

#### 2.1. Protocol:

22 male subjects , 11 with different knee abnormalities previously diagnosed by a professional. They undergo three movements to analyze the behavior associated with the knee muscle , gait , leg extension from a sitting position , and flexion of the leg up. The acquisition process was conducted with 4 electrodes ( Vastus Medialis , semitendinosus , biceps femoris and rectus femoris ) and the goniometer in the knee .

## 2.2. Instrumentation

Datalog equipment was used MWX8 by Biometrics of 8 digital channels and 4 analog channels , of which 4 for sampling were used SEMG and 1 for goniometry, these data were acquired directly to the computer MWX8 internal storage with microSD card and transmitted in Real-time Datalog software through bluetooth adapter , 14-bit resolution and sampling frequency of 1000Hz .

#### 2.3. Data configuration:

The total number of electrodes is 4, corresponding to the time series one for each channel (1 to 4). Each series contains  $\sim$  5 shares or motion repetitions for each subject.

#### 3. Attributes of information

Each data file contains 5 columns, organized as follows.

Segment		Lower Limb					
Channel	Ch1	Ch2	Ch3	Ch4	Ch5		
Muscle	RF	BF	VM	ST	FX		
Column	0	1	2	3	4		

**Segment:** defines the body part where the data are taken. **Channel:** corresponds to the electrode attached to a muscle.

Muscle: corresponds to the muscle being measured.

RF: Recto Femoral.
BF: Femoral Biceps.
VM: Vastus Medialis.
ST: Semitendinosus.
FX: Flexion at the knee.

# **Headings:**

File Name: 2Nsen.log

Channel 1: 'RF', 15300 values, engineering units: mV, no filters. Channel 2: 'BF', 15300 values, engineering units: mV, no filters. Channel 3: 'VM', 15300 values, engineering units: mV, no filters. Channel 4: 'ST', 15300 values, engineering units: mV, no filters.

Channel 5: 'FX', 765 values, engineering units: deg, no filters, extrapolated from 50

to 1000 samples per second.to 1000 samples per second.

#### 7. Number of classes:

The database contains 22 samples and 11 normal 11 with knee pathology, each subject has 3 different shots, one sitting, one standing and one gait.

## 8. Otros:

4 folders: A\_LOG, A\_TXT, N\_LOG Y N\_TXT.

Folder \_log contains data in .log format that can be loaded and analyzed by datalog software by Biometrics.

Folder \_txt contains SEMG data in columns and their headers